



Implant Compatibility Chart TITANIUM

Atlantis™ patient-specific abutments in Titanium and Atlantis GoldHue® are compatible with the implant systems and sizes listed below. For the most current implant compatibility, please visit www.atlantiscomp.com or contact Customer Service at (877)828-5628.

Implant Manufacturer	Implant Type and Diameter	Abutment Screw Interface Type
<i>Internal Connection Interfaces</i>		
Astra Tech	OsseoSpeed™ 4.0 ¹ , 3.5/4.0, 4.5/5.0	0.050" Hex Drive
Bio-Horizons®	3.5, 4.0, 5.0, 6.0	0.050" Hex Drive
Biomet 3i™	Certain® Microminiplant™ 3.4/3.25 Certain® 3.75/4.0, 5.0, 6.0, 4/5, 5/6	0.048" Hex Drive
Nobel Biocare®	Replace® NP (Narrow Platform 3.5) Replace RP (Regular Platform 4.3) Replace WP (Wide Platform 5.0, 6.0)	Nobel Biocare Unigrip®
Straumann™	Standard/Standard Plus Implant RN (Regular Neck) 4.8	Straumann Drive Compatible
Zimmer Dental®	Tapered Screw-Vent® 3.5, 4.5, 5.7 Screw-Vent® 3.5, 4.5 SwissPlus® 4.8	0.050" Hex Drive or 1.25mm Hex Drive
<i>External Hex Interfaces</i>		
Bio-Horizons®	3.5, 4.0, 5.0	0.050" Hex Drive
Biomet 3i™	Microminiplant™ 3.4/3.25 Miniplant™ 4.1/3.25 Standard 3.75, 4.0 Wide 5.0, 6.0 XP 3/4, 4/5, 5/6	0.048" Hex Drive
Innova Endopore® Innova Entegra™	3.5, 4.1, 5.0	0.050" Hex Drive
Lifecore®	Small Diameter 3.3, 3.4 Regular Diameter 3.75, 4.0 Wide Diameter 5.0, 5.5, 6.0	0.050" Hex Drive
Nobel Biocare®	Branemark® NP (Narrow Platform 3.3) Branemark® RP (Regular Platform 3.75, 4.0) Branemark® WP (Wide Platform 5.0, 5.5) SteriOss Replace® 3.1, 3.5, 4.3, 5.0, 6.0 SteriOss HL™ 3.25, 3.8, 4.5, 5.0, 6.0	Nobel Biocare Unigrip®
Sterngold-Impla-Med®	3.3, 3.75, 4.0, 5.0	0.050" Hex Drive

All orders for Atlantis Abutments are shipped with an Atlantis HT™ Titanium Abutment Screw at no extra charge. These abutment screws meet or exceed the implant manufacturer's recommendation for torque. **The implant manufacturer's recommended torque should always be used when placing the abutment screw.** Please consult the Atlantis *Instructions for Use* for more information.

¹Astra Tech "previous" 4.0 with packaging color-coded blue.



Implant Compatibility Chart

ZIRCONIA

Atlantis™ patient-specific abutments in Zirconia™ are compatible with the implant systems and sizes listed below. For the most current implant compatibility, please visit www.atlantiscmp.com or contact Customer Service at (877)828-5628.

Implant Manufacturer	Implant Type and Diameter	Abutment Screw Interface Type
<i>Internal Connection Interfaces</i>		
Astra Tech	OsseoSpeed™ 3.5/4.0, 4.5/5.0	0.050" Hex Drive
Biomet 3i™	Certain® Microminiplant™ 3.4/3.25 Certain® 3.75, 4.0, 5.0, 6.0, 4/5, 5/6	0.048" Hex Drive
Nobel Biocare®	Replace RP (Regular Platform 4.3) Replace WP (Wide Platform 5.0, 6.0)	Nobel Biocare Unigrip®
Zimmer Dental®	Tapered Screw-Vent® 3.5, 4.5, 5.7 Screw-Vent® 3.5, 4.5	0.050" Hex Drive or 1.25mm Hex Drive
<i>External Hex Interfaces</i>		
Bio-Horizons®	4.0, 5.0	0.050" Hex Drive
Biomet 3i™	Miniplant™ 4.1/3.25 Standard 3.75, 4.0 Wide 5.0, 6.0 XP 3/4, 4/5, 5/6	0.048" Hex Drive
Innova Endopore® Innova Entegra™	4.1, 5.0	0.050" Hex Drive
Lifecore®	Regular Diameter 3.75, 4.0 Wide Diameter 5.0, 5.5, 6.0	0.050" Hex Drive
Nobel Biocare®	Branemark® RP (Regular Platform 3.75, 4.0) Branemark® WP (Wide Platform 5.0, 5.5) SteriOss Replace® 4.3, 5.0, 6.0 SteriOss HL™ 3.8, 4.5, 5.0, 6.0	Nobel Biocare Unigrip®
Sterngold-Impla-Med®	3.75, 4.0, 5.0	0.050" Hex Drive

All orders for Atlantis Abutments in Zirconia are shipped with an Atlantis HTZ™ Titanium Abutment Screw at no extra charge. These abutment screws, which have been specifically designed for use with Atlantis Abutments in Zirconia, meet or exceed the implant manufacturer's recommendation for torque. **The implant manufacturer's recommended torque should always be used when placing the abutment screw.** Please consult the Atlantis *Instructions for Use* for more information.

All trademarks are the property of their respective owners.

This form was created by Astra Tech solely for use by its customers. Any other use is prohibited. ©2008 Astra Tech